

DOUBLE LUMEN OOCYTE ASPIRATION NEEDLE

ÇİFT LÜMENLİ YUMURTA TOPLAMA İĞNESİ

- ↳ 16G Double Lumen Needle maximize the quality of oocytes
- ↳ Echo marked tip gives extreme ultrasound definition
- ↳ Sharp needle tip ensures minimum trauma
- ↳ Silicone tap fits test tubes securely

Ref. No	Needle (G x cm)	Aspiration Line (cm)	Suction Line (cm)	Flushing Line (cm)
GOA1632	16 x 32	100	20	100
GOA1635	16 x 35	100	20	100

- ↳ 16G Çift Lümen İğne toplu alınan yumurta kalitesini artırır
- ↳ Ekolojik ucu ultrason altında mükemmel görünüm sağlar
- ↳ Keskin ucu travmaları minimuma indirir
- ↳ Silikon stoper, test tüplerine güvenli şekilde oturur