

STONE BASKET

TAŞ BASKET

- ↳ Stainless steel and nitinol alternatives
- ↳ Designed for endoscopic manipulation
- ↳ Detachable handle

- ↳ Paslanmaz çelik ve nitinol seçeneği
- ↳ Endoskopik kullanıma uygun dizayn
- ↳ Ayrılabilir sap



STAINLESS STEEL

Ref. No	Diameter (F)	Length (cm)	BASKET TYPE	WIRE TYPE
GSBH390	3	90	HELICAL	Round
GSBS390			PARALLEL	Flat
GSBH3110	3	110	HELICAL	Round
GSBS3110			PARALLEL	Flat
GSBH490	4	90	HELICAL	Round
GSBS490			PARALLEL	Flat
GSBH4110	4	110	HELICAL	Round
GSBS4110			PARALLEL	Flat
GSBH590	5	90	HELICAL	Round
GSBS590			PARALLEL	Flat
GSBH5110	5	110	HELICAL	Round
GSBS5110			PARALLEL	Flat

NITINOL

Ref. No	Diameter (F)	Length (cm)	BASKET TYPE	WIRE TYPE	Zero Tip Ref. No
GSBN1990	1.9	90	HELICAL	Round	GSBZ1990
GSBN19110	1.9	110	HELICAL	Round	GSBZ19110
GSBN19120	1.9	120	HELICAL	Round	GSBZ19120
GSBN2290	2.2	90	HELICAL	Round	GSBZ2290
GSBN22110	2.2	110	HELICAL	Round	GSBZ22110
GSBN390	3	90	HELICAL	Round	GSBZ390
GSBN3110	3	110	HELICAL	Round	GSBZ3110
GSBN490	4	90	HELICAL	Round	GSBZ490
GSBN4110	4	110	HELICAL	Round	GSBZ4110
GSBN590	5	90	HELICAL	Round	GSBZ590
GSBN5110	5	110	HELICAL	Round	GSBZ5110